


REFINED TECHNOLOGY

Premium Long Pulsed Laser System







QLEMENIC







MAKING A BETTER FUTURE

Based on innovative technology and professional workforce, We will strive to become one of the world's market leader rapidly growing.








info@elemenic.com






Laser Pioneering Company ELEMENIC



Premium Long pulsed Laser

Three wavelength Long pulsed laser in one machine – Alexandrite+Nd:YAG+KTP





Best-in-class specification

Big spot size (up to 24mm) & Power full energy (100J/300ms) & high speed(Long-width up to 3Hz, Short-width up to 10Hz)



Multi cooling system

Stronger and more effective with less side effect, also has button to choose from cryo cooling to air cooling





Halo-Toning


Halo-Toning uses special fractional hand piece and powerful long-pulse energy to reach the fine energy core by stimulating melanin and collagen with natural regenerative promotion effects which produces whitening and lifting effects at the same time.









Premium Long Pulsed Laser System - ELEMENT



Versatile Laser system & Multiple Treatment







"We believe quality is our most important value"

Making a Better Future

ELEMENIC researches and develops as well as manufactures products for laser business, medical and health care field. Based on innovative technology and professional workforce, We will strive to become one of the world's market leader rapidly growing.





SYSTEM SPECIFICATION



LONG PULSED LASER

Medium	Nd:YAG	Alexandrite	KTP
Wavelength	1064nm	755nm	532nm
Energy	100J	70J	10J
Pulse duration	0.1 - 300ms (multiple)		0.1 - 100ms
Pulse Rate	0.5 - 10Hz		
Spot size	2~24mm		2~10mm
Cooling system	Cryo Gas/Air cooler		
Beam Delivery	Lens-coupled optical fiber with Handpiece		
Dimensions	1000(H), 780(D), 400(W)		
Weight	130Kg		

PRODUCTS VIEW







SIDE VIEW



Office | 204, Singil-ro, Yeongdeungpo-gu, Seoul, Republic of Korea Factory | 3rd,34, Gyeongin-ro 3bun-gil, Bucheon-si, Gyeonggi-do, Republic of Korea Tel. 032 613-5007 Fax. 032 613 5009 E-mail. info@elemenic.com









